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- A. pH scouting for pre-concentration
- B. Immobilization of ligand on a chip
- C. Regeneration scouting
- D. Interaction analysis

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Getting started kit (anti-beta2u-globulin antibody, beta2u-globulin)

Product Name	Contents	Order code
Getting started kit	1mg/ml anti-beta2u-glubulin antibody 50 µl 100 µg/ml beta2u-glubulin) 50 µl	BR-1008-36
CM5 chip	Pack of 3ea Pack of 3ea	BR-1005-30 (For Biacore T200,T100) BR-1000-12 (For other system)
Amine coupling kit	750mg 1-ethy-3(3-dimethylaminopropyl)carbodiimide hydrochloride(EDC), 115mg N-Hydroxysuccinimide(NHS), 10.5ml 1.0M ethanolamine-HCl pH 8.5	BR-1000-50
Acetate 4.0	1 x50 ml	BR-1003-49
Acetate 4.5	1 x50 ml	BR-1003-50
Acetate 5.0	1 x50 ml	BR-1003-51
Acetate 5.5	1 x50 ml	BR-1003-52
Regeneration scouting kit	11 ml ethylene glycol 11 ml 10mM glycine-HCl pH1.5 11 ml 10mM glycine-HCl pH2.0 11 ml 10mM glycine-HCl pH2.5 11 ml 10mM glycine-HCl pH3.0 11 ml 4.0M magnesium chloride 11 ml 0.2M sodium hydroxide 11 ml 0.5% sodium dodecyl sulphate (SDS) 20 ml 10% Surfactant P20	BR-1005-56
HBS-EP+ 10X Buffer	1 x 1000ml 4 x 50ml	BR-1006-69 BR-1008-26

- Getting started kit (anti-beta2u-glubulin antibody, beta2u-glubulin)

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pH3.5 < Immobilization buffer < Ligand pI

3.1.2 Material

200mM HBS-P+ Buffer

10mM sodium acetate, pH 4.0, 4.5, 5.0, 5.5

50mM NaOH

Sensor Chip CM5

1mg/ml anti-beta2 μ -glubulin antibody (getting started kit)

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3.2.2 Material

Amino coupling kit (100mM NHS, 400mM EDC, 1M Ethanolamine)

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
Stoic ratio: if each analyte molecule can bind 5 ligand molecules, then the ratio is 1/5.

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<General guidance>


Ideal regeneration: Analyte response is consistent after repeated injections and within 10% of the level in the 1st injection

Too mild conditions: The analyte response decreases and the baseline response increases

Too harsh conditions: The analyte response decreases and the baseline response is constant or decreases

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
4. Reference

4.1 Using Biacore to Measure the Binding Kinetics of an Antibody-Antigen Interaction

Michael Murphy¹, Laure Jason-Moller¹, JoAnne Bruno¹

Biacore, Inc., Piscataway, New Jersey

Publication Name: Current Protocols in Protein Science

Unit Number: Unit 19.14

4.2 BIACORE Sensor Surface Handbook