


GE Healthcare
Life Sciences

Amersham™ CyDye™ fluors

Selection Guide




Welcome to Amersham CyDye fluors


Amersham CyDye fluors are versatile fluorophores designed for a broad range of applications. CyDye fluors combine exceptional photostability with bright and intense signals, resulting in high sensitivity and a highly reproducible performance at a biologically suitable pH range of 3 to 10. The use of multicolor CyDye fluors enables multiplexing without crosstalk. In addition, the CyDye fluors readily dissolve in aqueous buffers, which eliminates the need of organic solvents. Finally, the high purity and thorough QC testing of Amersham CyDye fluors generate high levels of chromophore and reactive dye content.

Narrow excitation and emission bands result in discrete signals from each fluor, which together with minimal cross-talk contributes to high accuracy.

CyDye emission profiles




Excitation and emission profiles of Cy3 and Cy5




Spectral properties of CyDye fluors

| | Color | Excitation max | Emission max |
|-------|--------|----------------|--------------|
| Cy2 | Green | 489 | 506 |
| Cy3 | Red | 550 | 570 |
| Cy3.5 | Red | 581 | 596 |
| Cy5 | Blue | 649 | 670 |
| Cy5.5 | Blue | 675 | 694 |
| Cy7 | Purple | 743 | 767 |




Multi-purpose labeling

- Wide selection of CyDye fluors for protein labeling for use in e.g., ELISA & immunoprecipitation
- Labeling strategy based on protein characteristics
 - Antibodies
 - Glycoproteins and carbohydrates
 - Amine-rich proteins




Fluorescent Western blotting

- Amersham ECL Plex™ for multiplexed, quantitative Western blotting
- Normalization to housekeeping protein in expression analysis
- Analysis of post-translational modifications and isoforms




2-D DIGE

- Size and charge matched CyDye DIGE fluors enabling multiplexed quantitative 2D electrophoresis
- Labeling strategy based on sample availability
- Saturation dyes
 - For scarce sample material
- Minimal dyes
 - When sample material is not limiting



Cell imaging

- CyDye labeled secondary antibodies for use in the IN Cell Analyzer Systems
 - Development of fixed cell assays for novel targets
 - Medium to high-throughput fluorescence imaging
 - Ready-to-use labeled antibodies enable multiplexed assays
- Streptavidin-conjugated CyDye fluors for use in
 - Flow-cytometry
 - In situ* hybridization



DNA, RNA, and oligo labeling

- For multicolor analysis in microarrays, aCGH, FISH, whole chromosome painting, karyotyping, and gene mapping
- A broad range of labeled nucleotides and reactive dyes available
- A range of CyDye labeled amidites are available for incorporation in oligosynthesis
- Provides options for direct labeling or post-labeling applications
- CyScribe™ labeling kits offer ease of use
- Custom pack sizes available upon request

CyDye fluor for a wide range of applications

Multi-purpose protein labeling

| | Quantity | Code no. | Cy2 | Cy3 | Cy3.5 | Cy5 | Cy5.5 | Cy7 |
|--|----------|----------|-----|-----|-------|-----|-------|-----|
|--|----------|----------|-----|-----|-------|-----|-------|-----|

5-packs

Mono-Reactive CyDye packs

Targeted and specific labeling of amine residues on antibodies, peptides, or oligonucleotides

| | | | | | | | | |
|------------|----------|---------|--|---|---|---|---|--|
| Cy3 mono | 5 × 1 mg | PA23001 | | • | | | | |
| Cy3.5 mono | 5 × 1 mg | PA23501 | | | • | | | |
| Cy5 mono | 5 × 1 mg | PA25001 | | | | • | | |
| Cy5.5 mono | 5 × 1 mg | PA25501 | | | | | • | |

Bis-Reactive CyDye packs


Suitable for general labeling of amine residues

| | | | | | | | | |
|-----------|----------|---------|---|---|---|---|---|--|
| Cy2 bis | 5 × 1 mg | PA22000 | • | | | | | |
| Cy3 bis | 5 × 1 mg | PA23000 | | • | | | | |
| Cy3.5 bis | 5 × 1 mg | PA23500 | | | • | | | |
| Cy5 bis | 5 × 1 mg | PA25000 | | | | • | | |
| Cy5.5 bis | 5 × 1 mg | PA25500 | | | | | • | |

Mono-Reactive Maleimide CyDye packs

A more selective method of antibody labeling that reduces the possibility of the label interfering with subsequent antibody-antigen reactions

| | | | | | | | | |
|--------------------------------------|----------|---------|--|---|--|---|--|--|
| Cy3 Maleimide mono-Reactive Dye Pack | 5 × 1 mg | PA23031 | | • | | | | |
| Cy5 Maleimide mono-Reactive Dye Pack | 5 × 1 mg | PA25031 | | | | • | | |



| | Quantity | Code no. | Cy2 | Cy3 | Cy3.5 | Cy5 | Cy5.5 | Cy7 |
|--|----------|----------|-----|-----|-------|-----|-------|-----|
|--|----------|----------|-----|-----|-------|-----|-------|-----|

Value packs

Mono-Reactive NHS ester

Targeted and specific labeling of amine residues on antibodies, peptides, or oligonucleotides

| | | | | | | | | |
|----------------------|-------|---------|--|---|---|---|---|---|
| Cy3 mono NHS ester | 1 mg | PA13101 | | • | | | | |
| | 5 mg | PA13105 | | • | | | | |
| | 10 mg | PA13104 | | • | | | | |
| | 25 mg | PA13106 | | • | | | | |
| | 50 mg | PA13102 | | • | | | | |
| Cy3.5 mono NHS ester | 1 mg | PA13601 | | | • | | | |
| | 5 mg | PA13605 | | | • | | | |
| | 25 mg | PA13606 | | | • | | | |
| | 50 mg | PA13602 | | | • | | | |
| Cy5 mono NHS ester | 1 mg | PA15101 | | | | • | | |
| | 5 mg | PA15100 | | | | • | | |
| | 10 mg | PA15104 | | | | • | | |
| | 25 mg | PA15106 | | | | • | | |
| | 50 mg | PA15102 | | | | • | | |
| Cy5.5 mono NHS ester | 1 mg | PA15601 | | | | | • | |
| | 5 mg | PA15605 | | | | | • | |
| | 10 mg | PA15604 | | | | | • | |
| | 25 mg | PA15606 | | | | | • | |
| | 50 mg | PA15602 | | | | | • | |
| Cy7 mono NHS ester | 1 mg | PA17101 | | | | | | • |
| | 5 mg | PA17105 | | | | | | • |
| | 10 mg | PA17104 | | | | | | • |
| | 50 mg | PA17102 | | | | | | • |

Bis-Reactive NHS ester

Suitable for general labeling of amine residues

| | | | | | | | | |
|---------------------|-------|---------|---|---|--|---|---|---|
| Cy2 Bis NHS ester | 5 mg | PA12000 | • | | | | | |
| | 10 mg | PA12004 | • | | | | | |
| | 50 mg | PA12002 | • | | | | | |
| Cy3 Bis NHS ester | 5 mg | PA13000 | | • | | | | |
| | 10 mg | PA13004 | | • | | | | |
| | 50 mg | PA13002 | | • | | | | |
| Cy5 Bis NHS ester | 5 mg | PA15000 | | | | • | | |
| | 10 mg | PA15004 | | | | • | | |
| | 50 mg | PA15002 | | | | • | | |
| Cy5.5 Bis NHS ester | 5 mg | PA15500 | | | | | • | |
| | 10 mg | PA15504 | | | | | • | |
| | 50 mg | PA15502 | | | | | • | |
| Cy7 Bis NHS ester | 5 mg | PA17000 | | | | | | • |
| | 10 mg | PA17004 | | | | | | • |
| | 50 mg | PA17002 | | | | | | • |

| | Quantity | Code no. | Cy2 | Cy3 | Cy3.5 | Cy5 | Cy5.5 | Cy7 |
|--|----------|----------|-----|-----|-------|-----|-------|-----|
|--|----------|----------|-----|-----|-------|-----|-------|-----|

Mono-Reactive Maleimide

A more selective method of antibody labeling that reduces the possibility of the label interfering with subsequent antibody-antigen reactions

| | | | | | | | | |
|----------------------|-------|---------|--|---|---|---|---|--|
| Cy3 mono maleimide | 1 mg | PA13131 | | • | | | | |
| | 5 mg | PA13130 | | • | | | | |
| | 25 mg | PA13136 | | • | | | | |
| Cy3.5 mono maleimide | 1 mg | PA13631 | | | • | | | |
| | 5 mg | PA13630 | | | • | | | |
| Cy5 mono maleimide | 1 mg | PA15131 | | | | • | | |
| | 5 mg | PA15130 | | | | • | | |
| | 25 mg | PA15136 | | | | • | | |
| Cy5.5 mono maleimide | 1 mg | PA15631 | | | | | • | |
| | 5 mg | PA15630 | | | | | • | |
| | 25 mg | PA15636 | | | | | • | |

Mono-Reactive Hydrazide

For the labeling of free carbonyl groups of glycoproteins and carbohydrates

| | | | | | | | | |
|----------------------|-------|---------|--|---|---|---|---|--|
| Cy3 mono hydrazide | 1 mg | PA13121 | | • | | | | |
| | 5 mg | PA13120 | | • | | | | |
| Cy3.5 mono hydrazide | 1 mg | PA13621 | | | • | | | |
| | 5 mg | PA13620 | | | • | | | |
| | 25 mg | PA13626 | | | • | | | |
| Cy5 mono hydrazide | 1 mg | PA15121 | | | | • | | |
| | 5 mg | PA15120 | | | | • | | |
| | 25 mg | PA15126 | | | | • | | |
| Cy5.5 mono hydrazide | 1 mg | PA15621 | | | | | • | |
| | 5 mg | PA15620 | | | | | • | |

Select-a-dye

Enables a variety of labeling strategies using NHS ester, maleimide, or hydrazide to be investigated

| | | | | | | | | |
|-----------------------|------------|---------|--|---|--|---|--|--|
| Cy3 Select-a-Dye pack | 3 × 0.5 mg | PA13123 | | • | | | | |
| Cy5 Select-a-Dye pack | 3 × 0.5 mg | PA15123 | | | | • | | |

Antibody labeling kits

A convenient way to tag antibodies with CyDye fluorescent dyes. Antibody kits contain all reagents required for 2 × 1 mg reactions

| | | | | | | | | |
|---------------------|--|---------|---|---|--|---|--|--|
| Cy2 Ab Labeling Kit | | PA32000 | • | | | | | |
| Cy3 Ab Labeling Kit | | PA33000 | | • | | | | |
| Cy5 Ab Labeling Kit | | PA35000 | | | | • | | |

Monoclonal antibody labeling kits

Monoclonal antibody kits contain all reagents required for 2 × 0.1 mg reactions

| | | | | | | | | |
|----------------------|--|---------|---|---|--|---|--|--|
| Cy2 mAb Labeling Kit | | PA32001 | • | | | | | |
| Cy3 mAb Labeling Kit | | PA33001 | | • | | | | |
| Cy5 mAb Labeling Kit | | PA35001 | | | | • | | |



Fluorescent Western Blotting

| Product | Code no. | Quantity | Cy2 | Cy3 | Cy5 |
|-----------------------------------------------------------------------------------|------------|----------|-----|-----|-----|
| ECL Plex goat-@-mouse IgG-Cy5, sufficient for 1000 cm ² membrane area | PA45009 | 150 µg | | | • |
| ECL Plex goat-@-mouse IgG-Cy5, sufficient for 4000 cm ² membrane area | PA45010 | 600 µg | | | • |
| ECL Plex goat-@-rabbit IgG-Cy5, sufficient for 1000 cm ² membrane area | PA45011 | 150 µg | | | • |
| ECL Plex goat-@-rabbit IgG-Cy5, sufficient for 4000 cm ² membrane area | PA45012 | 600 µg | | | • |
| ECL Plex goat-@-mouse IgG-Cy3, sufficient for 1000 cm ² membrane area | PA43009 | 150 µg | | • | |
| ECL Plex goat-@-mouse IgG-Cy3, sufficient for 4000 cm ² membrane area | PA43010 | 600 µg | | • | |
| ECL Plex goat-@-rabbit IgG-Cy3, sufficient for 1000 cm ² membrane area | 28-9011-06 | 150 µg | | • | |
| ECL Plex goat-@-rabbit IgG-Cy3, sufficient for 4000 cm ² membrane area | 28-9011-07 | 600 µg | | • | |
| ECL Plex goat-@-mouse IgG-Cy2, sufficient for 1000 cm ² membrane area | 28-9011-08 | 150 µg | • | | |
| ECL Plex goat-@-mouse IgG-Cy2, sufficient for 4000 cm ² membrane area | 28-9011-09 | 600 µg | • | | |
| ECL Plex goat-@-rabbit IgG-Cy2, sufficient for 1000 cm ² membrane area | 28-9011-10 | 150 µg | • | | |
| ECL Plex goat-@-rabbit IgG-Cy2, sufficient for 4000 cm ² membrane area | 28-9011-11 | 600 µg | • | | |



2-D DIGE

| | Code no. | Quantity | Minimal dye | Scarce sample |
|-----------------------------------------------------------------|------------|----------|-------------|---------------|
| CyDye DIGE Fluor, Cy2 minimal dye | 25-8010-82 | 5 nmol | • | |
| CyDye DIGE Fluor, Cy3 minimal dye | 25-8010-83 | 5 nmol | • | |
| CyDye DIGE Fluor, Cy5 minimal dye | 25-8010-85 | 5 nmol | • | |
| CyDye DIGE Fluor, Cy2 minimal dye | 25-8008-60 | 10 nmol | • | |
| CyDye DIGE Fluor, Cy3 minimal dye | 25-8008-61 | 10 nmol | • | |
| CyDye DIGE Fluor, Cy5 minimal dye | 25-8008-62 | 10 nmol | • | |
| CyDye DIGE Fluor, Cy2 minimal dye | RPK0272 | 25 nmol | • | |
| CyDye DIGE Fluor, Cy3 minimal dye | RPK0273 | 25 nmol | • | |
| CyDye DIGE Fluor, Cy5 minimal dye | RPK0275 | 25 nmol | • | |
| Minimal dye labeling kit | 28-9345-30 | 2 nmol | • | |
| CyDye DIGE Fluor, minimal labeling kit | 25-8010-65 | 5 nmol | • | |
| DIGE labeling kit for scarce samples | 25-8009-83 | 1 kit | | • |
| Scarce samples and preparative gel labeling | 25-8009-84 | 1 kit | | • |
| Preparative gel labeling for scarce samples, Cy3 saturation dye | 28-9366-83 | 400 nmol | | • |



Cell imaging

| | Code no. | Cy2 | Cy3 | Cy3.5 | Cy5 | Cy5.5 | Cy7 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-----|-----|-------|-----|-------|-----|
| CyDye labeled secondary antibodies | | | | | | | |
| <i>High content cell-based assays are powerful tools in drug discovery and cellular screening. Fluorescently labeled second antibodies allow researchers to develop their own fixed-cell assays for novel targets.</i> | | | | | | | |
| Cy2 labeled goat anti-mouse IgG, 1 mg | PA42002 | • | | | | | |
| Cy2 labeled goat anti-rabbit IgG, 1 mg | PA42004 | • | | | | | |
| Cy3 labeled goat anti-mouse IgG, 1 mg | PA43002 | | • | | | | |
| Cy3 labeled goat anti-rabbit IgG, 1 mg | PA43004 | | • | | | | |
| Cy5 labeled goat anti-mouse IgG, 1 mg | PA45002 | | | | • | | |
| Cy5 labeled goat anti-rabbit IgG, 1 mg | PA45004 | | | | • | | |

| Code no. | Cy2 | Cy3 | Cy3.5 | Cy5 | Cy5.5 | Cy7 |
|----------|-----|-----|-------|-----|-------|-----|
|----------|-----|-----|-------|-----|-------|-----|

Streptavidin-Fluor conjugates

Streptavidin conjugated to a range of fluorescent dyes for use in situ hybridization and flow cytometry

| | | | | | | |
|------------------|---------|---|---|---|--|--|
| Cy2-Streptavidin | PA42001 | • | | | | |
| Cy3-Streptavidin | PA43001 | | • | | | |
| Cy5-Streptavidin | PA45001 | | | • | | |

DNA and RNA labeling

| Feature | CyScribe™ First-Strand cDNA Labeling Kit | CyScribe Post-Labeling Kit | CyScribe aCGH Labeling kit |
|---------------------------------------------------------|------------------------------------------|----------------------------|---------------------------------|
| Signal brightness | ++ | +++ | +++ |
| Even incorporation of Cy3 and Cy5 | ++ | +++ | +++ |
| Starting material | mRNA or total RNA | mRNA or total RNA | DNA from cells or amplified DNA |
| Quantity of starting material | 5 to 25 µg of total RNA | 5 to 25 µg of total RNA | 1 µg of DNA |
| Quantity of starting material using mRNA | 50 ng to 1 µg | 100 to 500 ng | N/A |
| Possible to prepare a batch of unlabeled cDNA and store | No | Yes | No |
| Simplicity of protocol | +++ | + | ++++ |
| Time from DNA / RNA to labeled target | 3 h | 5.5 h | 3 h |
| Suitable for less experienced users | Yes | Less so | Yes |

+ Suitable
 ++ Good
 +++ Very Good
 ++++ Excellent

| Quantity | Code no. | Product information |
|----------|----------|---------------------|
|----------|----------|---------------------|

CyScribe labeling kits

For the generation of Cy3- and Cy5-labeled cDNA for microarray hybridization analysis

| | | | |
|--------------------------------------------------------------------------------------|--------------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| CyScribe First-Strand cDNA Labeling Kit | 25 reactions | RPN6200 | Provides reagents for preparation of Cy3 and Cy5 cDNA probes in first-strand cDNA synthesis reactions. Can be used with either Cy3 or Cy5 labeled dCTP and dUTP |
| CyScribe First-Strand cDNA Labeling System – dUTP | 50 reactions | RPN6201 | Contains x2 RPN6200, 25 nmol Cy3-dUTP and 25nmol Cy5-dUTP |
| CyScribe First-Strand cDNA Labeling System – dCTP | 50 reactions | RPN6202 | Contains x2 RPN6200, 25 nmol Cy3-dCTP and 25nmol Cy5-dCTP |
| CyScribe First-Strand cDNA Labeling System – dUTP with CyScribe GFX Purification Kit | 50 reactions | RPN6201X | Contains x2 RPN6200, 25 nmol Cy3-dUTP, 25nmol Cy5-dUTP and 50 Cyscribe GFX purification reagents |
| CyScribe First-Strand cDNA Labeling System – dCTP with CyScribe GFX Purification Kit | 50 reactions | RPN6202X | Contains x2 RPN6200, 25 nmol Cy3-dCTP, 25nmol Cy5-dCTP and 50 Cyscribe GFX purification reagents |
| CyScribe Post-Labeling Kit | 24 reactions | RPN5660 | Provides reagents for preparation of post-labeled Cy3 and Cy5 cDNA probes |
| CyScribe Array CGH Genomic DNA Labeling Kit | 30 reactions | 28909726 | For <i>in vitro</i> labeling of genomic DNA using Cy3 and Cy5 dCTP |
| CyScribe Array CGH Genomic DNA Labeling System | 30 reactions | 28919956 | Contains x1 28909726, plus x1 illustra GFX PCR DNA and Gel Band Purification Kit (28-9034-70) |

| | Quantity | Code no. | Product information |
|----------------------------------------------------------------------------------------------------------------|------------------------------------------------|------------|----------------------------------------------------------------|
| CyDye Fluorescent Nucleotides | | | |
| <i>CyDye labeled nucleotides are available in a range of bright, intense colors with narrow emission bands</i> | | | |
| Cy3-dCTP | 25 nmol | PA53021 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH7.0 |
| Cy3-dCTP | 25 nmol | PA53031 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy5-dCTP | 25 nmol | PA55021 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy5-dCTP | 250 nmol | PA55031 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy3-dUTP | 25 nmol | PA53022 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy3-dUTP | 250 nmol | PA53032 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy5-dUTP | 25 nmol | PA55022 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy5-dUTP | 250 nmol | PA55032 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy3.5-dCTP | 25 nmol | PA53521 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy5.5-dCTP | 250 nmol | PA55521 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH 7.0 |
| Cy3-CTP | 100 nmol | 25-8010-86 | Supplied as a 5 mM solution |
| Cy5-CTP | 100 nmol | 25-8010-87 | Supplied as a 5 mM solution |
| Cy3-UTP | 100 nmol | PA53026 | Supplied as a 5 mM solution |
| Cy5-UTP | 100 nmol | PA55026 | Supplied as a 5 mM solution |
| Value Pack dCTP | 5 × 25 nmol Cy3 dCTP + 5 × 25 nmol Cy5 dCTP | PA55321 | Bulk pack of labeled Cy3 and Cy5 labeled dCTP |
| Value Pack dUTP | 5 × 25 nmol Cy3 dUTP + 5 × 25 nmol Cy5 dUTP | PA55322 | Bulk pack of labeled Cy3 and Cy5 labeled dUTP |
| CyDye Post-Labeling Reactive Dye Pack | 12 × 40 000 pmol Cy3 + 12 × 40 000 pmol Cy5 | RPN5661 | Cy3 and Cy5 reactive NHS esters |
| Cy3 Post-Labeling Reactive Dye Pack | 12 × 40 000 pmol Cy3 | 25-8010-79 | Cy3 NHS esters |
| Cy5 Post-Labeling Reactive Dye Pack | 12 × 40 000 pmol Cy5 | 25-8010-80 | Cy5 NHS esters |

HyPer5™ Microarray Dye

| | | | |
|------------------------------------------------------------------------|---------------------------------------------------|------------|---------------------------------------------------------------|
| HyPer5 dCTP | 25 nmol | 28-9231-83 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH7.0 |
| HyPer5 dCTP | 250 nmol | 28-9231-84 | Supplied as a 1 mM solution, in 10 mM phosphate buffer, pH7.0 |
| Multipack containing 5 × 25 nmol Cy3 dCTP + 5 × 25 nmol HyPer5 dCTP | 5 × 25 nmol Cy3 dCTP + 5 × 25 nmol HyPer5 dCTP | 28-9231-85 | Cy3 and HyPer5 labeling dCTP |
| Cy3 and HyPer5 Post-Labeling Reactive Dye Pack | 12 × 40 000 pmol Cy3 + 12 × 15 000 pmol HyPer5 | 28-9224-19 | Cy3 and HyPer5 reactive NHS Esters |

Cy™ Amidites

For CyDye incorporation at any position in an oligonucleotide during synthesis

| | | | |
|---------------|--------|------------|--|
| Cy3 Amidite | 100 mg | 28-9172-98 | |
| Cy3 Amidite | 1 g | 28-9172-99 | |
| Cy5 Amidite | 100 mg | 28-9042-49 | |
| Cy5 Amidite | 1 g | 28-9021-58 | |
| Cy5.5 Amidite | 100 mg | 28-9042-50 | |
| Cy5.5 Amidite | 1 g | 28-9021-60 | |

Custom pack sizes are available for all Cy dye products. Please enquire with your local GE Sales representative for additional information.



Systems for fluorescence imaging

| | Cy2 | Cy3 | Cy5 | Comments |
|-------------------------------|-----|-----|-----|---------------------------------------------------------------|
| Single Fluorescence | | | | |
| ImageQuant™ LAS 4000 mini | • | | | |
| ImageQuant LAS 4000 | • | • | • | |
| ImageQuant LAS 4010 | • | • | • | |
| ImageQuant LAS 7000 | • | • | • | |
| Typhoon™ FLA 9500 | • | • | • | |
| Multiplex Fluorescence | | | | |
| ImageQuant LAS 4010 | (•) | • | • | Multiplex can be run as duplex either with Cy2/Cy5 or Cy3/Cy5 |
| Typhoon FLA 9500 | • | • | • | |
| DIGE | | | | |
| Typhoon FLA 9500 | • | • | • | |

For more information or to place your order, please visit www.gelifsciences.com or contact your local GE distributor.



For contact information for your local office,
please visit, www.gelifesciences.com/contact

www.gelifesciences.com/protein-purification

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden



GE, imagination at work and GE monogram are trademarks of General Electric Company.

Amersham, Cy, CyDye, CyScribe ECL Plex, GFX, HyPer5, ImageQuant, and Typhoon and are trademarks of GE Healthcare companies.

Cy3-dUTP or Cy5-dUTP, Cy3.5-dCTP or Cy5.5-dCTP, Cy3-dCTP or Cy5-dCTP: These products are manufactured for GE Healthcare UK Limited by Perkin Elmer Life Sciences under US patent numbers 5047519 and 5151507.

2-D Fluorescence Difference Gel Electrophoresis (2-D DIGE) technology is covered by US patent numbers 6,043,025, 6,127,134 and 6,426,190 and equivalent patents and patent applications in other countries and exclusively licensed from Carnegie Mellon University.

CyDye: This product is manufactured under an exclusive license from Carnegie Mellon University and is covered by US patent numbers 5,569,587 and 5,627,027.

Cy2: This product is manufactured under an exclusive license from Carnegie Mellon University and is covered by US patent numbers 5,486,616 and 5,627,027.

The purchase of CyDye products includes a limited license to use the CyDye products for internal research and development but not for any commercial purposes. A license to use the CyDye products for commercial purposes is subject to a separate license agreement with GE Healthcare. Commercial use shall include:

1. Sale, lease, license or other transfer of the material or any material derived or produced from it.
2. Sale, lease, license or other grant of rights to use this material or any material derived or produced from it.
3. Use of this material to perform services for a fee for third parties, including contract research and drug screening.

If you require a commercial license to use this material and do not have one, return this material unopened to GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden and any money paid for the material will be refunded.

© 2012 General Electric Company—All rights reserved.
First published Feb. 2012.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire, HP7 9NA
UK

GE Healthcare Europe, GmbH
Munzinger Strasse 5
D-79111 Freiburg
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327
Piscataway, NJ 08855-1327
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg., 3-25-1, Hyakunincho
Shinjuku-ku, Tokyo 169-0073
Japan