Biacore Maintenance Tools

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Table of Contents

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A. System Maintenance

1. Desorb	5
2. Desorb and Sanitize	6
3. Washing Buffer Tubing	9
4. Empty Buffer Tubing	11
5. Superclean	13
6. System Check	14
7. Flow System Wash	16

B. Chip Maintenance: Storing & Re-using



1. Wet storage	17
2. Dry storage	18

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< Maintenance 주기>

Interval time	Maintenance Tools
Daily & After Each Run	① Empty the waste bottle ② Desorb
Weekly	① Tube fitting and Pump leakage ② Desorb
Monthly	① Desorb and Sanitize ② Superclean ③ System Check
As Needed	① Washing Buffer Tubing ② Empty Buffer Tubing ③ Flow System Wash

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This procedure removes adsorbed material from the flow system.

Total run time is about 20 minutes.

Do not run this procedure below 20°C.

Note; Use the maintenance Chip or a used chip for this procedure. The ligand on the sensor chip may be damaged by the solutions used.

next


Required solutions(from Maintenance Kit);

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This procedure removes adsorbed material and disinfects the flow system. The procedure is divided into five steps. Total run time is about one hour followed by a recommended standby time of 3-4 hours.

Do not run this procedure below 20°C.

Do not abort this procedure after it is started.

Required solutions (from Maintenance kit)

BIAdesorb solution 1: one volume of 25ml and one volume of 15ml

BIAdesorb solution 2: one volume of 25ml and one volume of 15ml

BIAdisinfectant solution 1: one volume of 50ml and one volume of 30ml

Step 1

Place 25ml BIAdesorb solution 1 on the left hand tray and insert all four pump inlet tubes.

Place 15ml BIAdesorb solution 1 on right hand tray and insert the water inlet

Running Desorb and Sanitize procedure step 1. Please wait.

Time left: 00:08:50

Step 2

Place 25ml BIA desorb solution 2 on the left hand tray and insert all four pump inlet tubes.

Place 15ml BIA desorb solution 2 on right hand tray and insert the water inlet

Running Desorb and Sanitize procedure step 2. Please wait.

Time left: 00:08:50

Step 3

Wipe the pump inlet tubes with a moist tissue.

Place 50ml BIA disinfectant solution on the left hand tray and insert all four pump inlet tubes.

Place 30ml BIA disinfectant solution on the on right hand tray and insert the water inlet tube.

Running Desorb and Sanitize procedure step 3. Please wait.

Time left: 00:30:20

Step 4

Wipe the pump inlet tubes with a moist tissue.

Place water on the left hand tray and insert all four pump inlet tubes.

Place water on the on right hand tray and insert the water inlet tube.

Running Desorb and Sanitize procedure step 4. Please wait.

Time left: 00:03:30

Step 5

Place tube A in a HEPES or Tris buffer

Recommended concentration 10-50mmol/l.

Let tubes B, C and D hang in the air

Running Desorb and Sanitize procedure step 5. Please wait.

Time left: 00:13:30

The Desorb and Sanitize procedure is completed.

All the system to run in standby mode for at least 3-4 hours before performing a run

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This procedure removes adsorbed material from the flow system. The procedure is divided into steps. Total run time is about 30minutes.

NOTE: Used the maintenance Chi or a used chip for this procedure. The ligand on the sensor chip may be damaged by the solutions used.

Do not run this procedure below 20°C.
Do not abort this procedure after it is started.

Required solutions (from Maintenance Kit):

BIAdesorb solution 1. About 20ml

BIAdesorb solution 2. About 20ml

Step1

Place 20ml BIAdesorb solution 1 on the left hand tray and insert tube A.

Running wash buffer Tubing procedure step 1. Please wait.

Time left : 00:07:00

Step2

Wipe the tube with a moist tissue.

Place 20ml BIAdesorb solution 2 on the left hand tray and insert tube A

Running wash buffer Tubing procedure step 2. Please wait.

Time left : 00:07:00

Step3

Wipe the tube with a moist tissue.

Place buffer or water on the left hand tray and insert tube.

Running wash buffer Tubing procedure step 3. Please wait.

Time left : 00:12:20

The Wash Buffer Tubing procedure is completed.

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This procedure empties all four buffer selector inlet tubes. The procedure is divided into three steps. Total run time is about 20minutes.

Required solutions:
Deionized water
70% ethanol

Step 1

Place a bottle containing deionized water on the left hand plate and insert the four buffer inlet tubes.

Running Empty Buffer Tubing procedure step 1. Please wait.

Time left: 00:05:30

Step 2

Place a bottle containing at least 10ml 70% EtOH on the left hand plate and insert the four buffer inlet tubes.

Running Empty Buffer Tubing procedure step 2. Please wait.

Time left: 00:06:30

Step 3

Remove the tubes from the EtOH bottle and allow them to hang in the air.

This step empties the buffer selector of liquid.

Running Empty Buffer Tubing procedure step 3. Please wait.

Time left: 00:06:00

The Empty Buffer Tubing procedure is completed.

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This procedure cleans the flow system. Total run time is about 80 minutes.

Put tube A in filtered deionized warm water(40~50°C). 150ml

NOTE: Use the Maintenance Chip or a used chip for this procedure. The ligand on the sensor chip may be damaged by the solutions used.

Required solutions:

1% acetic acid

0.2M sodium bicarbonate

6M guanidine-HCl(or 50% DMSO)

10 mM HCl (or 10% DMSO)

If small molecule assays are run, exchange the last two wash solutions for DMSO solutions.

Running Superclean procedure. Please wait.

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Select test(s) to run.

This procedure should be run at 25°C with a new sensor Chip CM5 and with HBS-N as running buffer. Cheese Close if you need to change the sensor chip, reset the temperature or change running buffer.

If any injection is delayed, adjustments can be made by the software. The required test is injections and Refractometer.

Reagent pumps and blank injection Injections and Refractometer
 Mix Noise
 Merged and Dual injectios(optional) Buffer selector (optional)

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If the buffer selector shall be tested, put tube B, C and D in water.

If all inlet tubes are not intended to be used after System Check, don't forget to run Empty Buffer Tubing after the test.

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Select test(s)

A: Washing procedure for clogged IFC

B: Wash with Buffer

Wash solution : 50% Isopropanol/ 50mM NaOH/ 1M NaCl 1.2ml in r2a1. The chip may be affected. The run time is about 12minutes

This procedure washes the flow system with buffer at a high flow rate to clear any obstructions. Total run time is about 4 minutes

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1. Wet Maintenance

1.1. 재료: HBS-EP buffer (GE Healthcare) & A 50 ml polypropylene test tube

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